

MATERIAL TESTING SYSTEM

ACST series

Dynamic Compression Tester for Cushioning Materials

ACST-200

ACST can conduct the dynamic compression test to identify the characteristics of cushioning materials.

Shock response acceleration and displacement waveform of cushioning material are captured with ease when the compressive plate is dropped using this tester.

Features

- Easy operation with the automatic controller.
- Compressive plates for each mass are variable from 2.0kg to 50kg when in full spec.
- Velocity of drop plate also can be measured via optional indicator.
- With the specialized software, the characteristics of cushioning material is analyzed as stress-strain curve, cushion curve, cushion factor, etc. for practical use in the protective packaging design.

Compliance Standard


- JIS Z0235-2002
- ASTM D 1596-02 (2009)
- MIL-C-26861-87

Specification


Model	ACST-200
Specimen max size (mm)	220 × 220
Fall weight (kg)	Lightweight type : 2.3 ~ 10.0 Heavyweight type : 5.0 ~ 50.0
Max acceleration (m/s ²)	3000 (300G)
Max drop height (mm)	1200
Size (W × D × H mm)	W565 × D740 × H2730
Capacity (kg)	480
Controller Size (W × D × H mm)	W575 × D750 × H1220
Power supply	AC200V, 3-phases, 5A

Options


Velocity measuring device VS-2	Measurement of velocity of dropped weight
Analysis software for cushion characteristics CCA-Win	<ul style="list-style-type: none"> • Acceleration / Displacement - Time History • Cushion factor - Dynamic Stress Diagram • Maximum Acceleration - Static Stress Diagram • Maximum Displacement - Static Stress Diagram • Making of Test Data Record Table etc.



Acceleration and displacement wave form



Cushion Curve from some experimental data



Stress-Strain curve



SHINYEI TECHNOLOGY CO., LTD.

Tsukuba Office B47-11, Katoridai, Tsukuba, Ibaraki, 300-2657, JAPAN

Tel: +81-29-848-3571 FAX: +81-29-848-3572

International Sales Office

Shinyei BLDG. 5F, 77-1 Kyomachi, Chuo-ku, Kobe 651-0178, JAPAN

Tel: +81-78-392-6903 FAX: +81-78-332-1619